


Niti-S Esophageal Covered Stent [Beta-2]

for malignant and/or benign stricture and tracheoesophageal fistula



Features

- **Both distal and proximal retrieval strings help for easy removal or repositioning**
- **Outer double layers:** Silicone covered double layers prevent the risk of migration
- **Flexible and conformable structure:** Unfixed cell construction provides flexibility and conformability to fit in tortuous anatomy
- **Radiopaque marker:** 4 (four) at both ends & 3 (three) in the middle of each ring

Ordering Information

Distal Release				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
EK1808FNT2	18	8	20	70
EK1810FNT2		10		
EK1812FNT2		12		
EK1814FNT2		14		
EK1815FNT2		15		
EK1816FNT2		16		
EK1818FNT2		18		
EK1820FNT2		20		
EK2008FNT2	20	8	20	70
EK2010FNT2		10		
EK2012FNT2		12		
EK2014FNT2		14		
EK2015FNT2		15		
EK2016FNT2		16		
EK2018FNT2		18		
EK2020FNT2		20		
EK2208FNT2	22	8	20	70
EK2210FNT2		10		
EK2212FNT2		12		
EK2214FNT2		14		
EK2215FNT2		15		
EK2216FNT2		16		
EK2218FNT2		18		
EK2220FNT2		20		

Distal Release				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
EK2410FNT2	24	10	22	70
EK2412FNT2		12		
EK2414FNT2		14		
EK2415FNT2		15		
EK2416FNT2		16		
EK2418FNT2		18		
EK2420FNT2		20		
EK2610FNT2		26		
EK2612FNT2	12			
EK2614FNT2	14			
EK2615FNT2	15			
EK2616FNT2	16			
EK2618FNT2	18			
EK2620FNT2	20			
EK2810FNT2	28		10	22
EK2812FNT2		12		
EK2814FNT2		14		
EK2815FNT2		15		
EK2816FNT2		16		
EK2818FNT2		18		
EK2820FNT2	20			